

## OPTICLAR STAINLESS STEEL SPECULA TIDY

OPT700.000.030



### PRODUCT DESCRIPTION

## Stainless steel holder

Stainless steel wall mounted Tip Tidy

\*Handles, heads and specula not included\*

#### Includes:

- Stainless steel wall mounted Tip Tidy
- Holds 6 varieties of tips - metal or single patient use
- Easily cleaned and autoclaveable if required
- Bottom lip catches any residue from disinfected tips
- 2 large ports for holding instrument handles

### ADDITIONAL INFORMATION

<b>Animal</b>	N/A
<b>SKU</b>	OPT700.000.030
<b>Brand</b>	Opticlar

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights