

VETERINARY LED SLIT OTOSCOPE HEAD

OPT700.000.040



PRODUCT DESCRIPTION

Veterinary Slit Otoscope Head

Veterinary Slit Otoscope Head with Twist Lock Fitting.

Opticlar's VisionVet range is designed to simplify and enhance diagnostic examinations across species, from exotics to companion animals and larger breeds.

Engineered for durability and performance, each unit features TrueTone LED modules, delivering exceptional illumination with a long lifespan. The 10-year LED module guarantee means maintenance costs are significantly reduced over the lifespan of the equipment. The high-quality optics ensure excellent visualisation, while the robust construction offers lasting reliability.

Key Features

- Veterinary Slit Otoscope Head with Twist Lock Fitting
- Fits Opticlar and Welch Allyn handles as standard

Compatible with the following plastic and metal specula

Plastic

- OPTKS2.5V
- OPTKS3.5V
- OPTKS4.5V
- OPTKS5.5V
- OPT700.000.021
- OPT700.000.022
- OPT700.000.023
- OPT700.000.024
- OPT700.000.025

Metal

- OPT700.000.000-35
- OPT700.000.000-45
- OPT700.000.000-55
- OPT700.000.000-5L
- OPT700.000.001
- OPT700.000.001RG
- OPT700.000.002
- OPT700.000.002RG
- OPT700.000.003
- OPT700.000.004
- OPT700.000.007-2
- OPT700.000.010
- OPT700.000.011

Black

- OPT700.001MB
- OPT700.002MB
- OPT700.003MB

ADDITIONAL INFORMATION

Animal	Mixed
SKU	OPT700.000.040
Brand	Opticlar

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights