

OPTICLAR BATTERY POCKET DIAGNOSTIC SET - WITH MINI OPHTHALMOSCOPE AND SLIT OTOSCOPE

OPT700.010.010D



PRODUCT DESCRIPTION

Opticlar P2 Mini Pro Veterinary Diagnostic Set - Plastic Tips - 100 Hour Run Time

Opticlar's VisionVet range is designed to simplify and enhance diagnostic examinations across species, from exotics to companion animals and larger breeds.

Engineered for portability, durability and performance, each unit features TrueTone LED modules, delivering exceptional illumination with a long lifespan. The 10-year LED module guarantee means maintenance costs are significantly reduced over the lifespan of the equipment. The high-quality optics ensure excellent visualisation, while the robust construction offers lasting reliability.

An automatic on/off function and compatibility with traditional manufacturer connections help protect your existing equipment investments. Additionally, the multi-voltage design allows seamless switching between cell power and rechargeable options without the need to replace bulbs.

Set Includes:

- ProVet LED slit otoscope
- Pocket Pro ophthalmoscope
- 2 x P2 polymer pocket handles
- 5 x plastic tips - 57.0 x 4.0mm diameter
- 5 x plastic tips - 65.0mm x 6.0mm diameter
- Stylish pouch

Key Features

- The basic student diagnostic set features the ProVet LED slit head otoscope and Pocket Pro ophthalmoscope
- Ideal for a second domiciliary set, to leave in vehicles or as a backup surgery set

- Excellent LED illumination
- Handles manufactured from impact-resistant polymer
- Integral on-off switch
- Supplied with 10 x single-use plastic tips

P2 Mini Pro Pocket Handle - AA Battery

- Takes X2 AA batteries
- Acrow Thread
- Running Time: 100 hours of continuous use
- Rotary control allows the LED brightness to be easily adjusted

Pocket Slit Otoscope compatible with the following plastic and metal specula

Plastic

- OPTKS2.5V
- OPTKS3.5V
- OPTKS4.5V
- OPTKS5.5V
- OPT700.000.021
- OPT700.000.022
- OPT700.000.023
- OPT700.000.024
- OPT700.000.025

Metal

- OPT700.000.000-35
- OPT700.000.000-45
- OPT700.000.000-55
- OPT700.000.000-5L
- OPT700.000.001
- OPT700.000.001RG
- OPT700.000.002
- OPT700.000.002RG
- OPT700.000.003
- OPT700.000.004
- OPT700.000.007-2
- OPT700.000.010
- OPT700.000.011

Black

- OPT700.001MB
- OPT700.002MB
- OPT700.003MB

ADDITIONAL INFORMATION

Animal	Mixed
SKU	OPT700.010.010D
Brand	Opticlar
Diagnostic Head Type	Co-axial Ophthalmoscope, Slit Otoscope

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights