

OPTICLAR FINOFF TRANSILLUMINATOR

OPT100.000.013



PRODUCT DESCRIPTION

Veterinary Finoff transilluminator

Opticlar's Finoff transilluminator illuminates the fundus via the sclera, and is often referred to as the 'muscle light'. Designed for use when there is a presence of intraocular tumours in patients where cloudy media is suspected.

Key Features

- LED illumination source
- Finoff Transilluminator, for scleral transillumination and checking pupil reflex
- Suitable for frontal and maxillary sinus examination
- Fits any Opticlar or Welch Allyn handle with a twist lock connection
- Fits to Keeler and Heine handles by using an adaptor
- Ideal for transillumination and for use with an Indirect Ophthalmoscopy lens
- Suitable for fluorescence examinations.

*Please note: Handle not included

ADDITIONAL INFORMATION

| | |
|---------------|----------------|
| Animal | Mixed |
| SKU | OPT100.000.013 |
| Brand | Opticlar |

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights