

OPTICLAR BI-VALVE SPECULUM

OPT100.000.015



PRODUCT DESCRIPTION

Veterinary Bi-Valve Spreadable speculum

Opticlar's VisionVet range is designed to simplify and enhance diagnostic examinations across species, from exotics to companion animals and larger breeds.

Engineered for durability and performance during nasal, rectal and vaginal procedures.

Key Features

- Ideal for various examinations including nasal, rectal and vaginal
- Fits any Opticlar or Welch Allyn handle with twist lock connection
- Heine or Keeler handle compatible with adapter
- Fingertip control to facilitate the easy removal of foreign bodies LED illumination and fibre optic light stem

ADDITIONAL INFORMATION

Animal	Mixed
SKU	OPT100.000.015
Brand	Opticlar

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights