

OPTICLAR OPERATING OTOSCOPE HEAD ONLY

OPT100.000.012



PRODUCT DESCRIPTION

Veterinary Operating Otoscope - Head Only.

Key Features

- True Tone LED illumination source
- Open design for easier instrument access. Swivelling lens for easier instrumentation.
- Excellent LED illumination
- 3.0x magnification from scratch proof window lens for enhanced viewing
- Distal fibre optics gives cool light to examination site with no reflections and provides excellent view of the tympanic membrane
- Swivel viewing lens cannot be lost and is guarded on lower edge from potential breakage
- Wide angle viewing window for easy instrumentation under magnification
- Head available with German or American tip fitting
- OpenHeadDesign - allows introduction of instruments without obscuring the view

Compatible with the following plastic specula:

- WA22002
- WA22003
- WA22004
- WA22005
- WA22009
- WA22064
- WA22065
- WA22067
- WA22069
- WA22100
- WA22160

ADDITIONAL INFORMATION

Animal

Mixed

SKU

OPT100.000.012

Brand

Opticlar

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights