

OPTICLAR HANDLE - LARYNGOSCOPE LITHIUM ION

OPT600.090.010

**PRODUCT DESCRIPTION****Opticlar Rechargeable handle with LED illumination.**

A hygienic handle provides a comfortable grip yet is easier to clean than a standard diagnostic handle.

Key Features

- Rechargeable Lithium-Ion Battery
- The whole handle is autoclavable
- Fits Miller or Macintosh blades
- High-performance LED light
- Recharge handle via Opticlar docking station or mini USB cable

ADDITIONAL INFORMATION

Animal	Mixed
SKU	OPT600.090.010
Brand	Opticlar

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights