

OPTICLAR METAL SPECULA

OPT700.000.001-GROUPED



PRODUCT DESCRIPTION

Opticlar Veterinary Metal Specula - Small / Medium / Large

For use with Veterinary Slit & Pocket Otoscopes

- Small, 57.0 x 4.0mm (OPT700.000.001)
- Medium, 57.0 x 4.0mm (OPT700.000.002)
- Large, 57.0 x 4.0mm (OPT700.000.003)

Standard Specula - Supplied as standard in the Veterinary slit otoscope set and available individually. The bayonet fitting is compatible with Opticlar and Heine Slit Otoscope heads. Fully autoclaveable

Specifications:

- Metal construction
- Autoclaveable
- Bite-resistant

ADDITIONAL INFORMATION

SKU

OPT700.000.001-grouped

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights