

OPTICLAR STANDARD SPECULA - MATT BLACK

OPT700.000.001MB-GROUPED

**PRODUCT DESCRIPTION****Opticlar Standard Specula - MATT BLACK - Small / Medium / Large****For use with Veterinary Slit & Pocket Oscopes**

- Opticlar Standard Specula - Matt Black - Small - 57mm x 4.0mm (OPT700.000.001MB)
- Opticlar Standard Specula - Matt Black - Medium - 65mm x 6.0mm (OPT700.000.002MB)
- Opticlar Metal Specula - Matt Black - Large - 90mm x 7.0mm (OPT700.000.003MB)

Standard Specula - Matt Black - The Matt Black finish completely eliminates glare. Supplied as standard in the Veterinary slit otoscope set and available individually. The bayonet fitting is compatible with Opticlar and Heine Slit Oscope heads. Fully autoclaveable

Specifications:

- Metal construction
- Autoclaveable
- Bite-resistant

ADDITIONAL INFORMATION**SKU**

OPT700.000.001MB-grouped

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights