


LOCKABOX ONE™ MEDI

LOCKABOX



PRODUCT DESCRIPTION

The Lockabox One™ Medi

Designed with a versatile range of healthcare uses in mind for whoever – and whatever – you need to keep safe. Its tough, shatter-resistant polycarbonate and three-digit combination lock can securely store medicines and more.

- Made with food-grade and medical-grade materials
- High-impact and shatter-resistant polycarbonate construction
- Changeable 3-digit combination code lock
- Patented clamshell design
- Stackable and can be drilled and fixed as required
- Suitable for fridges and freezers (down to -80°C)
- Vents to allow airflow
- Opaque design
- Made in Britain from recyclable polymer
- Fully assembled and ready to use
- 2-year warranty

Please note: Colour may vary

Lockabox One™ Dimensions:

External: 31cm (L) x 21cm (W) x 17 cm (H) (including vents and lock mechanism)

Internal: 28.3 cm (L) x 19.2cm (W) x 16.1 cm (H)

Door Opening: 17.5cm (W) x 15.3cm (H)

ADDITIONAL INFORMATION

SKU	LOCKABOX
Lock type	Key Lock

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.
Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!
Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.
Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights