

MOBILE INFUSION STAND

I4110017



PRODUCT DESCRIPTION

Easy to assemble mobile infusion (drip) stand*

Part of the Burtons value range, this practice essential with its chrome finish looks great and is easy moved around the practice thanks to its lightweight construction and 5 star castor base.

- Height adjustable (135-215 cm)
- Manoeuvrable and lightweight
- Chrome finish
- 4 Hook design
- 5 star castor base

***Infusion Pump Not Included**

ADDITIONAL INFORMATION

| | |
|---------------------|----------|
| Animal | Mixed |
| SKU | I4110017 |
| Colour/Style | Silver |
| Brand | Burtons |

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights