

MOBILE ENDOSCOPY CART WITH 4 SHELVES

END-CTR900N



PRODUCT DESCRIPTION

A mobile 4-shelf endoscopy cart with 4-inch lockable footbrake wheels

Ideal for holding your complete endoscopy stack in one unit, the cart is made primarily from an aluminium and steel plate frame and has ABS/PC easy-clean shelf surfaces. It is complete with 2 adjustable endoscopy holders.

Key Features:

- 4 shelves
- 1 Arm
- 4-inch lockable castor wheels
- 2 x Adjustable Flexible Endoscopes Probe Brackets (Type C)

Technical Specification:

- Weight: 32kg
- Dimensions (HxWxD): 1183mm x 473mm x 592mm
- Max Load: 10kg
- Outer Shelf Size Dimensions (WxL): 436mm x 403mm
- Shelf Depth: 30mm
- Shelf Partition (WxH): 403mm x 240mm
- Material: Aluminium, Steel Plate, ABS/PC

ADDITIONAL INFORMATION

SKU	END-CTR900N
Animal	N/A
Overall Product Dimensions (mm)	1183H X 473W X 592D
Generic Size	Large

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.
Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!
Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.
Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights