

BURTONS EN-V5+ INFUSION PUMP & WARMING TUBE BUNDLE

EN-V5-TUBE-BUNDLE



PRODUCT DESCRIPTION

Burtons EN-V5+ Infusion Pump and Warming Tube Bundle

Introducing the Burton's EN-V5+ Infusion Pump, a groundbreaking device engineered to streamline veterinary care and enhance patient safety. Setting a new industry standard, the EN-V5+ is uniquely designed to store up to 200 multiple giving set calibrations, eliminating the need for separate pumps for different applications like paediatric, standard infusion, or blood transfusions. This unparalleled versatility not only saves costs but also simplifies your clinic's equipment management.

Experience unparalleled ease of use with the EN-V5+'s large, intuitive touchscreen. Effortlessly adjust settings and modes, and clearly monitor infusion progress, including a precise 'Time remaining' display for efficient patient management.

Designed for the dynamic veterinary environment, its powerful battery offers up to 10 hours of continuous operation, providing seamless portability between theatre and recovery. The included interchangeable IV pole bracket and cage mount, with a quick-release Allen key fitting, ensure secure and adaptable placement. The exclusive Burton's edition kit further enhances convenience with an 'L' shaped power lead, expertly designed to keep cords safely away from curious patients. Your EN-V5+ comes pre-calibrated to Burton's Giving Sets, with custom calibration available upon order.

The Burton's EN-V5 Infusion pump has been manufacturer-tested and is suitable for administering whole blood, plasma, and other blood products. However, gravity feeding remains the gold standard, as the use of pumps may carry a risk of premature haemolysis. Clinicians should apply their own judgement for each patient and follow appropriate guidance. Care should be taken to use correctly calibrated giving sets, recognising that some residual risk may remain.

Key Features for Superior Care

Multi-Calibration Storage: Store up to 200 giving set calibrations for ultimate versatility.

Comprehensive Kit: Includes interchangeable IV pole bracket, cage bracket, and right-angle power cord as standard.

- **Night Mode:** Dims screen brightness for reduced patient stress in dark recovery environments
- **Bloods Compatible:** Safely suitable for blood transfusions
- **Enhanced Safety:** Right-angle power lead protects against patient interference
- **Touchscreen Control:** Intuitive parameter settings and editing directly on the screen
- **Flexible Infusion:** 5 infusion modes plus micro mode control, with rates from 0.1ml/h to 2000ml/h
- **IP44 Splashproof:** Durable construction for demanding veterinary settings
- **Robust Alarms:** Comprehensive alarm system for VTBI, occlusion, pressure, battery, air bubbles, and more

Integrate Infusion Warming with Ease

Enhance patient comfort and maintain healthy body temperatures with the included EN-V5WT Warming Lead, exclusively compatible with the EN-V5+ pump. This warming lead conveniently holsters into the side of the unit and operates up to 40°C, providing improved coverage and heat distribution across the IV giving set. Note: Warming function is operational when running from mains power.

ADDITIONAL INFORMATION

Animal	Mixed
SKU	EN-V5-TUBE-BUNDLE
Category Sub Heading	Burtons EN-V5+ Infusion Pump with Warming Tube
Battery	Yes
Screen Type	Digital Display, Full Colour, Touchscreen
Volume Control	Yes
Flow Rate	0.1 - 2000 ml/h

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.
Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!
Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.
Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights