

BURTONS EN-V3 INFUSION PUMP

EN-V3



PRODUCT DESCRIPTION

The EN-V3 Infusion Pump is designed specifically for small to medium-sized animals, offering accurate IV delivery in a user-friendly, compact design. Its bright 2.8" colour screen and intuitive menu buttons ensure quick, easy operation for everyday veterinary use. Compatible with all Burton's giving sets and supporting multiple infusion modes, the EN-V3 provides straightforward, safe, and reliable patient infusion.

Key Features

- Multiple infusion modes: Offers ml/h (Rate), Time, Body Weight, Drip, and Micro modes for all treatment protocols
- High accuracy: Delivers infusion rates ranging from 0.1 to 200ml/h with an accuracy of $\pm 5\%$
- Adjustable KVO (Keep Vein Open): From 0.1 to 5ml/h, ensuring minimal fluid retention
- Occlusion detection: Four occlusion levels (150-900mmHg) to prevent overpressure and ensure patient safety, plus an occlusion alarm for when pressure exceeds preset levels
- Last therapy function: Quickly resumes previous infusion settings for rapid treatment initiation
- Advanced patient safety: The EN-V3 is equipped with multiple alarm systems to ensure patient safety, including low battery, system errors, VTBI and volume alarms for when the set volume is nearing completion.
- Brackets included: Pole and Cage brackets included
- Right-angled power cable: for keeping cables away from prying paws

Technical Specifications

- Flow Rate Range: 0.1-2000 ml/h (minimum increment: 0.01 ml/h)
- VTBI (Volume to Be Infused): 0.01-9999.99
- Bolus Rate: 0.1-2000 ml/h (both programmed and manual bolus available)
- Battery Life: Lithium Polymer battery (2600mAh) provides over 5 hours of operation at 25 ml/h;
- Dimensions: 138 mm (L) x 131.5 mm (W) x 90 mm (H)
- Weight: Approximately 1.5 kg (including pole clamp)
- Power Supply: AC 110V-240V, 50/60Hz; DC 10-16V

ADDITIONAL INFORMATION

Animal	Small Animal
SKU	EN-V3
Overall Product Dimensions (mm)	138mmL x 131.5mmW x 90mmH
Generic Size	Small
Colour/Style	Blue
Brand	Burtons
Battery	Yes
Screen Type	Full Colour
Volume Control	Yes
Flow Rate	0.1-2000 ml/h (minimum increment: 0.01 ml/h)

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights