

QUICKDROP 0.25 LEAD APRON (FEMALE)

146-856-F-M-L

**PRODUCT DESCRIPTION****Medium / Long 0.25 Lead Apron in Royal Blue with Velcro Adjustment (Female)**

Durable, lightweight lead aprons designed to protect veterinary professionals during radiographic procedures.

Combining effective radiation shielding with comfort and ease of movement, these aprons offer reliable protection without compromising on performance or practicality.

The protective material used in radiation aprons and accessories contains heavy metals that reduce your exposure to harmful scatter radiation, shielding veterinary professionals from the harmful effects of radiation exposure.

Key Features

- 2-Year warranty on apron
- Lightweight Lead Material
- With Pocket
- Velcro Strap Adjusters

Technical Specification (Medium / Long Apron Size)

- Chest Size: 94cm – 101cm (58cm W)
- Waist Size: 81cm
- Hip Size: 101cm
- Length: 106cm (Long)

ADDITIONAL INFORMATION

Animal	N/A
SKU	146-856-F-M-L
Overall Product Dimensions (mm)	Medium/Long
Brand	Burtons

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.
Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!
Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.
Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights