

## DRAMINSKI STANDOFF FOR LINEAR PROBE

DRAM-BLUE-STAND-OFF



### PRODUCT DESCRIPTION

## The standoff allows for the examination of tendons using an abdominal linear probe

For use with the Draminski high-frequency linear probe.

### ADDITIONAL INFORMATION

|  |                      |
|--|----------------------|
| <b>Animal</b>                          | N/A                  |
| <b>SKU</b>                             | DRAM-BLUE-STAND-OFF  |
| <b>Brand</b>                           | Draminski            |
| <b>ultrasound probes compatibility</b> | Draminski BLUE       |
| <b>Ultrasound Probe Uses</b>           | Large Animal Tendons |

#### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

