

DRAMINSKI LOW FREQUENCY CONVEX PROBE

DRAM-CON-LOW



PRODUCT DESCRIPTION

For pregnancy diagnostics in sheep and goats. Abdomen and thorax of cattle and horses. Horse musculoskeletal system. Umbilical region in calves.

Specifications:

- Frequency range: **2.0 - 5.0 MHz**
- Central frequency: **3.5 MHz**
- Scanning depth: **60 - 350 mm**
- Lens radius: **60 mm**
- Number of crystals: **128**

ADDITIONAL INFORMATION

| | |
|--|---|
| Animal | Equine and large animal |
| SKU | DRAM-CON-LOW |
| ultrasound probes compatibility | Draminski BLUE |
| Ultrasound Probe Uses | Pregnancy diagnostics in sheep and goats. Abdomen and Thorax of cattle and horses. Horse musculoskeletal system. Umbilical region in calves |

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights