

## DRAMINSKI ABDOMINAL HIGH FREQUENCY LINEAR PROBE

DRAM-LINEAR-HIGH



### PRODUCT DESCRIPTION

---

**For tendon examination in a horse or lung and abdomen examination of small animals**

#### Specifications

- Frequency range: 6.0 - 14.0 MHz
- Central frequency: 10.0 MHz
- Scanning depth: 40 - 100 mm
- Probe face length: 40 mm
- Number of crystals: 128

### ADDITIONAL INFORMATION

<b>Animal</b>	Companion animals, Equine and Large Animal
<b>SKU</b>	DRAM-LINEAR-HIGH
<b>ultrasound probes compatibility</b>	Draminski BLUE
<b>Ultrasound Probe Uses</b>	Equine tendon examinations, Small animal lung and abdomen examinations

## Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights