

WHEELED STAND FOR TENDONS EXAMINATION

DRAM-WHEELED-STAND



PRODUCT DESCRIPTION

The tendon examination stand ensures the right position during examinations of large animal tendons and facilitates movement around the patient

- Easily adjust the height and viewing angle of the scanner.
- Equipped with one probe holder
- A quick-folding system makes it easy to transport in a vehicle.

Specifications

- The height of the stand measured at the center of the ultrasound display when the stand is tilted at 45°: - 45 cm
- The height of the stand measured at the center of the ultrasound display with the stand fully extended: - 57 cm
- Dimensions when folded (without ultrasound scanner): 38 cm x 40 cm x 18 cm
- Dimensions when folded (with ultrasound scanner): **41 cm x 40 cm x 24 cm**
- Weight (without ultrasound scanner): 4.5 kg
- Number of swivel wheels: 4 (all wheels with lock)

ADDITIONAL INFORMATION

SKU	DRAM-WHEELED-STAND
Brand	Draminski
Animal	Equine
Overall Product Dimensions (mm)	38 cm x 40 cm x 18 cm when folded without scanner

Generic Size

Large

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.
Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!
Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.
Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights