

## ETHICON ENSEAL HANDPIECES

ETH-NSL-BUNDLE\*



### PRODUCT DESCRIPTION

## The ENSEAL X1 Tissue Sealers are bipolar electrocautery instruments for use with the Ethicon GEN11 generator.

They are ideal for use during open or laparoscopic surgical procedures to cut and seal vessels, and to cut, grasp and dissect tissue during surgery.

These semi-disposable handpieces must be sterilised with ethylene oxide between uses. Our service team can arrange sterilisation for you - please call [01622 834350](tel:01622834350).

### ADDITIONAL INFORMATION

<b>Animal</b>	Mixed
<b>SKU</b>	ETH-NSL-BUNDLE*

#### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

