

## ETHICON GEN11 VESSEL SEALING AND CUTTING SYSTEM

ETH-GEN11-BUNDLE\*



### PRODUCT DESCRIPTION

**The Ethicon GEN11 Generator offers veterinary practices advanced yet intuitive vessel sealing and cutting capabilities.**

**The energy control system is designed to support both HARMONIC™ Ultrasonic and ENSEAL™ Advanced Bipolar devices.**

Perfect for laparoscopic procedures, the GEN11 is powered by Adaptive Tissue Technology that intelligently senses changes in tissue and device conditions to deliver the optimal amount of energy, enhancing precision and protecting surrounding tissue during every procedure.

It really is as simple as plugging in your selected handpiece, following the onscreen prompts and using the trigger handpiece to activate.

The large, high-resolution touchscreen ensures effortless setup and operation, with clear, easy-to-read visuals from multiple angles. Automatic instrument recognition allows seamless connection of compatible devices, reducing setup time and potential errors in the theatre.

Compact and designed with the modern operating environment in mind, the GEN11 fits easily into any theatre layout and will easily fit on your laparoscopy trolley. Software is upgradeable via USB, ensuring your investment stays ready for future innovations.

#### Key Features:

- Supports all current and next-generation HARMONIC® and ENSEAL® instruments
- Adaptive Tissue Technology for intelligent, efficient energy delivery
- The intelligent system makes it incredibly quick and easy to use
- Automatic device recognition
- High-resolution, broad-view, touchscreen interface
- Compact design for space-saving, fits onto an endoscopy trolley seamlessly
- Software updates available to extend functionality over time

At Burtons, the Ethicon GEN11 Generator is part of our commitment to providing veterinary teams with reliable, cutting-edge surgical technology backed by expert advice and service. To discover how the GEN11 can integrate seamlessly into your practice, speak with our experts today.

Please note: all semi-disposable handpieces must be sterilised with ethylene oxide between uses. Our service team can arrange sterilisation for you - please call 01622 834350

## ADDITIONAL INFORMATION

<b>Animal</b>	Mixed
<b>SKU</b>	ETH-GEN11-BUNDLE*

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.  
**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!  
**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.  
**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights