

## AIRFLOW A2 HEATED INSUFFLATOR

END-300-900-021



## PRODUCT DESCRIPTION

### Features:

- User-friendly interface
  - Highest safety standard for patients
  - Pre-insufflation mode for use with Veress Cannulas
  - High gasflow rates
- Delivery contains:
- Main Unit
  - Sterile Filter END-300-902-007
  - Wrench Set END-300-902-006
  - CO2 High Pressure Tube END-300-902-003
  - Standard Silicone Tube END-300-902-005
  - Mains Cable
  - Operating Manual

## ADDITIONAL INFORMATION

<b>Animal</b>	Companion Animals, Exotics
<b>SKU</b>	END-300-900-021
<b>Brand</b>	SPITZ
<b>Procedure Type</b>	Bronchoscopy, Cystoscopy, Respiratory Examination, Rhinoscopy

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights