

5.5MM SELF-CUTTING MULTIPORT TROCAR SYSTEM

END-300-226-055



PRODUCT DESCRIPTION

MultiPort Trocar System 5.5mm self-cutting, threaded

- Reusable and autoclavable seals
- Easy opening by pushing the lever
- Stainless Steel stopcock
- Ergonomic Stainless Steel housing
- Self-centering sleeves

ADDITIONAL INFORMATION

Animal	Companion Animals, Exotics
SKU	END-300-226-055
Brand	SPiTZ
Procedure Type	Bronchoscopy, Cystoscopy, Respiratory Examination, Rhinoscopy
Endoscopy Device	Accessory
Animal Type	Companion Animal

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights