

PATIENT NIBP HOSE

SUN-91-0028*BUNDLE



PRODUCT DESCRIPTION

NIBP hose specifically for use with SunTech patient monitors, compatible with the SunTech Vet 20/25/30 in lengths of 4ft, 6ft, and 10ft. Simply slot and twist into back of your SunTech machine, secure fitting, and easy to clean.

ADDITIONAL INFORMATION

SKU	SUN-91-0028*Bundle
Accessory Parameter	NIBP
Monitor Compatibility	SunTech Vet20, SunTech Vet20E, SunTech Vet25, SunTech Vet25E, SunTech Vet30, SunTech Vet30E
Change Frequency	When Broken

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights