

## LABCOLD PHARMACY BENCHTOP REFRIGERATOR - 66L

124-005-CONFIG\*



### PRODUCT DESCRIPTION

---

## A pharmacy refrigerator featuring the latest Labcold Intellicold technology

Designed for the safe storage of cold chain medicines in line with Public Health England guidance. The 24 range features an enhanced touchscreen controller that delivers accurate timekeeping, intelligent event logging, and Bluetooth connectivity. Providing faster access to temperature data than ever before via the Labcold Connect Android App. No more having to manually upload temperature logs to your practice PC! This flat, easy-to-use controller automatically locks when not in use to prevent accidental changes.

Three independent internal temperature probes provide advanced temperature readings, ensuring that temperature readings accurately reflect the true product conditions throughout the fridge. This provides accurate, meaningful alarms when medicines are at risk. Reducing the risk of oversized or unrepresentative loads causing inaccurate readings. These all include a low-temperature safeguard that prevents freezing.

Each fridge is supplied with a FREE UK 5-year warranty, once registered, for complete peace of mind.

A UKAS Temperature calibration certificate gives you peace of mind that your temperatures are calibrated correctly. It can also be used if any health and safety issues arise around the correct storage of drugs, and can be shown to an auditor during inspections. These certificates can be purchased as an optional extra at the time of purchase.

#### Key features:

- Fitted with the IntelliCold touch screen controller
- 24/7 automatic temperature data logging via Bluetooth (Labcold Android App)
- Event-based logging
- Fan-assisted cooling
- Door lock options – key lock with two keys, or digital lock
- T3 Sensory & low temperature cut
- Realistic pre-fitted simulated load
- Controller and alarm battery back-up in the event of a power cut
- Minimum/maximum temperature recording and display

- Auto defrost
- Fan-assisted cooling
- Access port for additional monitoring equipment
- Environmentally friendly refrigerant
- Self-closing door
- Internal light
- Remote alarm contacts
- Brackets for wall mounting
- Compliant with all leading professional guidelines
- **Alarms:**
  - Door open
  - Power failure
  - High/ low temperature

**Specifications:**

- Capacity: 66L
- Style: Benchtop / wall mounted
- UK Warranty: 5 years
- External dimensions (mm): 735h x 450w x 510d
- Internal dimensions (mm): 588h x 331w x 226d
- Temperature range: 2°C to 8 °C
- Fridge set point: 5°C
- Internal material: Moulded plastic
- External material: Powder-coated steel
- Shelves: 3
- Weight: 34kg
- Door style: Solid
- Lock Type: Choice of key or digital lock
- Number of temperature probes: 3
- Power consumption: Kwh/24hrs

**ADDITIONAL INFORMATION**

|  |                               |
|--|-------------------------------|
| <b>SKU</b>                             | 124-005-config*               |
| <b>Overall Product Dimensions (mm)</b> | 735 x 450 x 510               |
| <b>Generic Size</b>                    | Small                         |
| <b>Brand</b>                           | Labcold                       |
| <b>Catagory Sub Heading</b>            | Choice of key or digital lock |
| <b>Alarm</b>                           | Yes                           |
| <b>Fridge Capacity</b>                 | 66L                           |
| <b>Door Type</b>                       | Regular                       |

## Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights