

AIRWAY ADAPTOR SIDESTREAM CAPNOGRAPHS - STRAIGHT FOR COMPATIBLE WITH

129-860-3



PRODUCT DESCRIPTION

Airway Adaptor Sidestream Capnographs - Straight . Change every other month. Compatible with M880B, Edan X12, Mindray uMEC, Mindray ePM12 & SurgiVet Advisor monitors.

ADDITIONAL INFORMATION

Animal	Mixed
SKU	129-860-3
Accessory Parameter	CO ²
Monitor Compatibility	BLT M880B, Edan X12, Mindray ePM12, Mindray uMEC, SurgiVet Advisor
Change Frequency	Every Other Month

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

