

**BP SPHYGMOMANOMETER**

129-900

**PRODUCT DESCRIPTION**

**Designed for ease and accuracy, this sphygmomanometer is ideal for quick and reliable blood pressure assessments in veterinary settings**

Its combined bellows and gauge tube design allows comfortable one-handed operation, while the clear 48mm dial ensures effortless reading of pressure values.

Featuring a straightforward 1-tube system with a push-fit cuff connection, this unit keeps setup simple and efficient.

**Key Features:**

- 1 tube system
- Combined bellows and gauge tube
- Can be used with any slip luer/ grip cuffs or Doppler
- 48mm dial for clear reading of pressure values
- Pre-set for right-handed use but can easily be switched to left-handed

**ADDITIONAL INFORMATION**

<b>Animal</b>	N/A
<b>SKU</b>	129-900
<b>Accessory Parameter</b>	NIBP
<b>Monitor Compatibility</b>	VetBP

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights