

EDAN IM50/60/70 SERIES - ECG CABLE,5 LEAD,SNAP,DEFIB,IEC,3.5M,REUSABLE

E-01.57.471471



PRODUCT DESCRIPTION

Edan IM50/60/70 Series - ECG Cable,5 lead,snap,Defib,IEC,3.5m,reusable. Change when broken.

ADDITIONAL INFORMATION

Animal	N/A
SKU	E-01.57.471471
Accessory Parameter	ECG
Change Frequency	When Broken

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.
Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!
Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.
Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights