

MINDRAY OESOPHAGEAL ECG & TEMP PROBE SUITABLE FOR SMALL ANIMALS UNDER 13KG

UMEC-TCN21-0037-00



PRODUCT DESCRIPTION

Mindray Oesophageal ECG & Temp Probe suitable for Small Animals under 13kg. Change when broken. Compatible with Mindray uMEC & Mindray ePM12 monitors.

ADDITIONAL INFORMATION

Animal	N/A
SKU	UMEC-TCN21-0037-00
Brand	Mindray
Accessory Parameter	ECG Temp
Monitor Compatibility	Mindray ePM12, Mindray uMEC
Change Frequency	When Broken

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.
Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!
Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.
Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

