

## REUSABLE MINDRAY COMPATIBLE SKIN SURFACE TEMP

IMEC-0011-30-37393



### PRODUCT DESCRIPTION

Reusable Mindray Compatible Skin Surface Temp. Change when broken. Compatible with Mindray uMEC & ePM12 monitors.

### ADDITIONAL INFORMATION

|                              |  |
|------------------------------|--|
| <b>Animal</b>                | N/A  |
| <b>SKU</b>                   | IMEC-0011-30-37393   |
| <b>Brand</b>                 | Mindray  |
| <b>Accessory Parameter</b>   | Temp   |
| <b>Monitor Compatibility</b> | BLT M880B, Edan X12, Mindray ePM12, Mindray uMEC, SurgiVet Advisor |
| <b>Change Frequency</b>      | When Broken  |

#### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

