

## SAMPLE TUBE FOR THE M880 CAPNOGRAPH

129-860-2

### PRODUCT DESCRIPTION

---

Sample Tube For the M880 Capnograph. Change evry 2-4 weeks. Compatible with M880B, Edan X12, Mindray uMEC, Mindray ePM12 & SurgiVet Advisor monitors.

### ADDITIONAL INFORMATION

<b>Animal</b>	N/A
<b>SKU</b>	129-860-2
<b>Accessory Parameter</b>	CO <sup>2</sup>
<b>Monitor Compatibility</b>	BLT M880B, Edan X12, Mindray ePM12, Mindray uMEC, SurgiVet Advisor
<b>Change Frequency</b>	2-4 weeks

#### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.  
**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!  
**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.  
**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights