

MACH LED 130 LAMP WITH CEILING MOUNTED KIT

MACHLED130-C



PRODUCT DESCRIPTION

Mach LED 130 Lamp with Ceiling Mounted Kit

Producing a cool light the LED technology is much more effective than conventional light sources. The temperature increase in the surgeons working area is almost non-existent. Features include a faceted multi-lens system, a colour rendering index of Ra=95 and superior handling design allowing the light to be positioned exactly to the wound field.

When the ceiling height is below 2.8m or above 3.0m please state the fixation height and height to underside of the false ceiling so the correct length pole can be supplied.

Specifications Common To All Models		
Technical Specifications	Mach LED 130 F	Mach LED 130
Light Intensity at 1 Meter Distance (lux)	60,000	50,000
Colour Rendering Index R(a) General	95	95
Light-Field Size (cm)	14 - 25 Focusable	Fixed Focus 17
Colour Temperature	4,500	4,500
Electronic Light Intensity Control On Lamp Head	50 - 100%	50 - 100%
Power Consumption	28 W	28 W
Number of LED's	19	19
Mains	230 V 50 Hz	230 V 50 Hz
Working Distance (cm)	70 - 140	70 - 140
Diameter of Light Head	33	33
Height Adjustment (cm)	123	123
Weight (Ceiling, Wall, Mobile)	15 Kg, 8 Kg, 20Kg	

ADDITIONAL INFORMATION

Animal	N/A
SKU	MACHLED130-C
Autoclavable Handle	No
Focusable	No
Head Type	Single
Light Source	LED
Lux @1 meter	60,000
Procedures	Dental, Exam, Minor-Operation, Task

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights