

SUNTECH ACCUVET 9.8FT (3M) TEMPERATURE PROBE

SUN-52-0023-00



PRODUCT DESCRIPTION

SunTech AccuVet 9.8ft (3m) Temperature Probe. Change when broken. SunTech Vet30 & SunTech Vet30E monitors.

ADDITIONAL INFORMATION

Animal	N/A
SKU	SUN-52-0023-00
Brand	SunTech
Accessory Parameter	Temp
Monitor Compatibility	SunTech Vet20, SunTech Vet20E, SunTech Vet25, SunTech Vet25E, SunTech Vet30, SunTech Vet30E
Change Frequency	When Broken

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.
Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!
Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.
Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

