

VETRONIC BLOOD PRESSURE EXTN LINE 200CM/1-2.0MM FOR LIGHTNING

VSL-71100.2



PRODUCT DESCRIPTION

Vetronic Blood Pressure Extn Line 200cm/1-2.0mm for Lightning. Single Use. Compatible with Vetronic Lightning monitors.

ADDITIONAL INFORMATION

Animal	N/A
SKU	VSL-71100.2
Brand	Vetronic
Accessory Parameter	IBP
Monitor Compatibility	Vetronic Lightning
Change Frequency	When Broken

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.
Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!
Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.
Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

