

ADAPTABLE ANAESTHETIC MACHINE SYSTEM BUNDLE

102-068-PROMO



PRODUCT DESCRIPTION

A classic from the Burtons Adaptable range

Bundle includes adaptable frame, Burtons vaporiser, vaporiser fitting kit, mini-lack, cyclo-flo, and concentrator.

Built upon the original Burtons AAA (Adaptable Anaesthesia Apparatus) frame concept, the Adaptable Anaesthesia System is designed for versatility and can form the basis of a variety of anaesthesia circuits.

Perfect for piped gas, this machine also accommodates two F-size oxygen cylinders in hinged clamps on the rear and an on-board oxygen concentrator in the base tray. Any and all of these oxygen supply methods can be used without the need to modify your machine.

For peace of mind in theatre, the Adaptable includes a number of carefully considered safety features. In cases of power outage, the manual changeover valve allows you to quickly and seamlessly switch from your on-board concentrator to your backup cylinders. The emergency flush button contains two integral safety valves for pressure relief and to prevent expired gases flowing back through the circuit to the vaporiser and beyond.

Key Features:

- Complete Burtons Lifetime stainless steel frame – durable, robust, easy-clean;
- High-accuracy 0-10l/min logarithmic rotameter;
- Emergency O₂ flush includes integral pressure relief and non-return safety valves;
- Connects seamlessly with the piped oxygen in your practice, F-size cylinders, or with an on-board concentrator;
- Rapid manual change-over valve for switching between concentrator and cylinder oxygen supplies – vital for power outages;
- Includes a bullnose pressure regulator for your F-size cylinder, with quick-release Schrader fittings;
- Can accommodate up to two vaporisers;
- Fully mobile, solid four-wheel base with independent wheel locks;
- Full warranty, supported by Burtons' nationwide Service Team;
- Designed and manufactured by Burtons in the UK;
- Optional height-adjustable monitor shelf also available; concentrator and vaporiser(s) sold separately.

Technical specs:

- Dimensions: 1380 x 460 x 570mm HWD
- Weight: c20kg

ADDITIONAL INFORMATION

Animal	Small Animal
SKU	102-068-promo
Overall Product Dimensions (mm)	43 x 37 x 65 WDH
Gas Type	Concentrator
Gas Supply	'F' Size Cylinder(s)
Flowmeter Type	Flowmeter
Concentrator	Yes
Optional O2 Alarm	No
Optional N2O	No
Emergency O2 Flush	No

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.
Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!
Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.
Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights