

BIOPSY TRANSFER CASSETTES

END-BTC



PRODUCT DESCRIPTION

Biopsy embedding cases for processing and embedding samples.

Mainly used for the handling, storage and transportation of biopsies and small tissues. Featuring fine mesh holes to ensure fluid pass through and proper drainage. As well as a 30-degree angle for easy labelling. They are resistant to solvents and chemical agents.

To seal the cassettes, the lid is broken off and locked in the lower section. They can be opened and securely sealed multiple times.

Key features:

- Solvent resistant
- 25 per pack
- Large labelling areas on the front and sides
- With lid attached
- Surface perforated with holes
- Dimensions: 40.5mmL x 28.5mm x 7mmD

ADDITIONAL INFORMATION

SKU	END-BTC
Animal	Mixed

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights