

BURTONS VETCLAVE 23B AUTOCLAVE (23L B-CLASS)

23B-VETCLAVE



PRODUCT DESCRIPTION

Reliable, High-Capacity Sterilisation for Veterinary Practices

The Burton's VetClave 23B Autoclave has been designed to meet the daily demands of busy practices, delivering consistent, dependable performance you can rely on. Combining an extra-large chamber with intuitive operation, it helps your team maintain the highest standards of decontamination[GK1.1] while keeping workflows running smoothly.

Engineered with one of the deepest chambers in its class, its 530 mm depth allows for the effective processing of larger instrument trays and longer surgical tools in a single cycle. This reduces the need for repeat loads and saves valuable time. Fast, optimised cycles mean sterile, dry instruments can be ready in as little as 31 minutes, ensuring you always have the equipment you need, exactly when you need it.

Designed with usability, safety and compliance in mind, the VetClave 23B simplifies the everyday sterilisation of hollow, porous and solid loads. With enhanced drying performance and antimicrobial surface protection, it helps ensure instruments are ready to use quickly while maintaining strong infection control standards throughout your practice.

Key Features

Extra-deep 23L chamber

One of the deepest chambers in its class (530 mm), providing the flexibility to accommodate larger instrument trays and longer surgical tools in a single cycle, ideal for busy practices.

Fast, optimised sterilisation cycles

Delivers sterile, dry instruments in as little as 31 minutes, helping to reduce downtime between procedures and maintain a consistent workflow throughout the day.

Integrated water management system

Designed for ease and convenience, the built-in frontal water pump enables quick, mess-free filling directly from ground-level containers at the touch of a button.

Enhanced Drying Performance

Customisable additional drying cycles (5 or 15 minutes) ensure instruments are completely moisture-free, ready for immediate use or safe storage.

Advanced Safety Features

Includes automatic door locking, continuous water quality monitoring and maintenance alerts to ensure safe, reliable operation and consistent sterilisation results.

Antimicrobial Surface Protection

BioCote® antimicrobial coating helps reduce microbes on surfaces by up to 99.99%, supporting a cleaner, more hygienic clinical environment.

Designed for Efficiency

The VetClave 23B is engineered to support the demands of modern veterinary practices, helping you optimise time, resources and workflow efficiency:

- Multiple sterilisation cycles tailored to different instrument types and loads
- Energy-saving standby mode and delayed start functionality
- Efficient water and power usage to reduce running costs
- Built for long-term reliability with a robust, durable design and Burtons equipment maintenance plans available

Cycle Suite

- 134°C Porous Load Wrapped
- 121°C Porous Load Wrapped
- 134°C Hollow Load Wrapped
- 121°C Hollow Load Wrapped
- 134°C Solid Load Unwrapped
- 121°C Solid Load Unwrapped
- Steam Penetration Vacuum Test – Leak Detection
- Helix / B&D Test

Specifications

- Chamber Size: 236 × 530 mm
- Capacity (B-Class):
 - o 3.5 kg wrapped instruments
 - o 4.5 kg unwrapped instruments
 - o 1.5 kg textiles
- Device Dimensions (W × D × H): 505 × 690 × 400 mm
- Power Supply: 220–240 V, 50–60 Hz
- Power Consumption: 2400 W Water Tank Capacity: 4 L (clean & used)
- Weight: 58 kg

What's Included

- 4 × trays
- Reversible rack and tray holder
- Fill and drain tubes
- Sponge

ADDITIONAL INFORMATION

SKU	23B-VETCLAVE
Overall Product Dimensions (mm)	505 × 690 × 400 mm
Generic Size	Medium
Colour/Style	Stainless Steel
Brand	Burtons
Number of Cycles	2
Vacuum	No
Chamber Size	22 Litres
Autoclave Accessory/Spares	Yes

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights