

MACH LED 110 EXAMINATION SPOT LIGHT

MACHLED110*BUNDLE



PRODUCT DESCRIPTION

Basic and affordable

The Mach LED 110 is small and compact, it is easily fitted to desks and ideal for additional lighting in examinations and consults

Key Features

- Flexible gooseneck arm - Easy to manoeuvre
- Ultra compact design, ideal for small spaces
- LED technology has a minimum lifespan of 40,000 hours
- Authentic colour rendering to reduce eyestrain
- 30,000 lux at 0.5 m distance
- Colour rendering Index R(a) 95
- Total power consumption: 7 watts
- Colour temperature (Kelvin): 4000

ADDITIONAL INFORMATION

Animal	Mixed
SKU	machled110*bundle
Colour/Style	White
Brand	Dr Mach
Autoclavable Handle	No
Battery Backup	No

Colour Control	No
Colour Temperature (k)	4000
Focusable	No
Head Type	Single
LED Life (hours)	40,000
Light Source	LED
Lux @1 meter	15,000 lux
Proceduredure	Exam, Task

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights