

## MACH LED 110 EXAMINATION SPOT LIGHT

MACHLED110\*BUNDLE



### PRODUCT DESCRIPTION

## Basic and affordable

The Mach LED 110 is small and compact, it is easily fitted to desks and ideal for additional lighting in examinations and consults

#### Key Features

- Flexible gooseneck arm - Easy to manoeuvre
- Ultra compact design, ideal for small spaces
- LED technology has a minimum lifespan of 40,000 hours
- Authentic colour rendering to reduce eyestrain
- 30,000 lux at 0.5 m distance
- Colour rendering Index R(a) 95
- Total power consumption: 7 watts
- Colour temperature (Kelvin): 4000

### ADDITIONAL INFORMATION

<b>Animal</b>	Mixed
<b>SKU</b>	machled110*bundle
<b>Colour/Style</b>	White
<b>Brand</b>	Dr Mach
<b>Autoclavable Handle</b>	No
<b>Battery Backup</b>	No

<b>Colour Control</b>	No
<b>Colour Temperature (k)</b>	4000
<b>Focusable</b>	No
<b>Head Type</b>	Single
<b>LED Life (hours)</b>	40,000
<b>Light Source</b>	LED
<b>Lux @1 meter</b>	15,000 lux
<b>Proceduredure</b>	Exam, Task

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights