

## MACH LED 115 EXAMINATION SPOT LIGHT

MACHLED115-BUNDLE



### PRODUCT DESCRIPTION

---

**The Mach Led 115 desk mounted spot light is sturdy and flexible with 60,000 lux at a distance of 0.5m. Making it ideal for minor procedures.**

#### Features

- Ultra compact design, ideal for small rooms.
- LED technology has a minimum life-span of 40,000 hours

- Easy to manoeuvre flexible gooseneck arm
- Authentic colour rendering to reduce eye strain (rendering index R(a) 95)
- Desk clamp option
- 60,000 lux at 0.5m distance
- Light head: 22cm
- Number of LEDs: 7
- Total power consumption: 10 watts

## ADDITIONAL INFORMATION

<b>Animal</b>	N/A
<b>SKU</b>	MACHLED115-bundle
<b>Colour/Style</b>	White
<b>Brand</b>	Dr Mach
<b>Autoclavable Handle</b>	No
<b>Battery Backup</b>	No
<b>Colour Control</b>	No
<b>Colour Temperature (k)</b>	4000
<b>Focusable</b>	No
<b>Head Type</b>	Single
<b>LED Life (hours)</b>	40,000
<b>Light Source</b>	LED
<b>Lux @1 meter</b>	30,000
<b>Procedure</b>	Exam, Task

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

