

MACH LED 115 EXAMINATION SPOT LIGHT

MACHLED115-BUNDLE



PRODUCT DESCRIPTION

The Mach Led 115 desk mounted spot light is sturdy and flexible with 60,000 lux at a distance of 0.5m. Making it ideal for minor procedures.

Features

- Ultra compact design, ideal for small rooms.
- LED technology has a minimum life-span of 40,000 hours

- Easy to manoeuvre flexible gooseneck arm
- Authentic colour rendering to reduce eye strain (rendering index R(a) 95)
- Desk clamp option
- 60,000 lux at 0.5m distance
- Light head: 22cm
- Number of LEDs: 7
- Total power consumption: 10 watts

ADDITIONAL INFORMATION

Animal	N/A
SKU	MACHLED115-bundle
Colour/Style	White
Brand	Dr Mach
Autoclavable Handle	No
Battery Backup	No
Colour Control	No
Colour Temperature (k)	4000
Focusable	No
Head Type	Single
LED Life (hours)	40,000
Light Source	LED
Lux @1 meter	30,000
Procedure	Exam, Task

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.
Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!
Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.
Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

