

MACH LED 130 PLUS EXAMINATION LIGHT

MACHLED130PLUS-CONFIG



PRODUCT DESCRIPTION

Ultra Powerful Mach LED 130 Plus Examination Light

The Mach LED 130 Plus is a compact examination light with powerful 100,000* lux LED technology. Producing a cool light the LED technology is much more effective than conventional light sources. The temperature increase in the working area is almost non-existent. Features include a faceted multi lens system, a colour rendering index of Ra=95 and superior handling design allowing the light to be positioned exactly to the wound field.

Key Features:

- 100,000 Lux (at a distance of 1m)
- Facetted multi lens system
- Colour rendering index of Ra=95
- Authentic colour rendering to reduce eye strain
- Superior handling design - allowing light to be positioned where needed
- 0.5 temperature increase of light head
- Fixed focus
- Colour temperature (Kelvin): 4300
- Dimming range (%): 50-100
- Number of LED: 19
- Size of the light field (cm): 12 (fix)
- Total power consumption (W): 28
- Life-span of the LED's (h) 50.000
- Working distance (in cm): 70 - 140
- Chassis of the light (cm): 33

Hygienic

High quality materials provide security and avoid infiltrations of dust, dirt and liquids inside the lamp head. The disk sealings of the light outlets and the circumferential sealing cord guarantee a perfect product down to the smallest detail.

Effective

Powerful LED-technology helps keeping a cool head through an almost inexistent heat development in the surgeon's head area and wound field.

Innovative

Facetted multi-lens system guarantees the lowest shadiness and absolute homogeneity. A multitude of computer-calculated facetted lenses create punctual illuminated light fields, increasing the contrast effect and therefore exact illumination of the wound field.

ADDITIONAL INFORMATION

Animal	N/A
SKU	MACHLED130PLUS-config
Colour/Style	White
Brand	Dr Mach
Autoclaveable Handle	No
Battery Backup	No
Colour Control	No
Colour Temperature (k)	4300
Head Type	Single
LED Life (hours)	50,000
Light Source	LED
Lux @1 meter	100,000

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights