

MACH LED 130F FOCUSABLE EXAMINATION LIGHT

MACHLED130F-CONFIG



PRODUCT DESCRIPTION

Part of the Mach LED 130 range, the Mach LED 130F is a compact, focusable examination light with adjustable light field sizes.

Producing a cool light the LED technology is much more effective than conventional light sources, with favourable flow channels that prevent disruptive heat development. The temperature increases in the surgeon's working area is almost non-existent. Features include a faceted multi lens system, a colour rendering index of Ra=96 and superior handling design allowing the light to be positioned exactly over the wound field. The precise light field control with the focusable controls via the handle help to concentrate focus and ascertain specific diagnostics during surgical tasks.

Available with optional focusable head allowing easy viewing of the working field and also with an optional autoclavable and easy to sterilise handle on the Mobile and Wall Mounted models.

When the ceiling height is below 2.8m or above 3.0m please state the fixation height and height to underside of the false ceiling so the correct length pole can be supplied.

Key Features:

- Brightness (LUX): 70,000 (at a distance of 1m)
- Faceted multi lens system for homogenous and shadow-free illumination
- Colour Rendering Index (CRI): Ra = 96 for accurate tissue colour identification
- Authentic colour rendering to reduce eye strain
- Superior handling design - allowing light to be positioned where needed
- 0.5°C temperature increase
- Colour temperature (Kelvin): 4500
- Dimming range (%): 50-100
- Number of LED: 19
- Size of the light field (cm): 13 - 19 (focusable)

- Total power consumption (W): 16
- Lifespan of the LED's (h): 60.000
- Working distance (in cm): 70-140
- Chassis of the light (cm): 33

Hygienic

High quality materials provide security and avoid infiltrations of dust, dirt and liquids inside the lamp head. The disk sealings of the light outlets and the circumferential sealing cord guarantee a perfect product down to the smallest detail.

Effective

Powerful LED-technology helps keeping a cool head through an almost non-existent heat development.

ADDITIONAL INFORMATION

Animal	N/A
SKU	MACHLED130F-config
Colour/Style	White
Brand	Dr Mach
Autoclavable Handle	Yes
Battery Backup	No
Colour Control	No
Colour Temperature (k)	4500
Focusable	Yes
Head Type	Single
LED Life (hours)	60,000
Light Source	LED
Lux @1 meter	70000
Produceure	Dental, Exam, Minor-Operation

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

