

## WELCH ALLYN PORTABLE ZIP UP CASE SETS

WA962-CASE-SET



### PRODUCT DESCRIPTION

## Welch Allyn Diagnostic Zip Up Case Set - Designed for veterinary professionals who do not have access to constant mains power.

3.5V Veterinary Diagnostic Set including coaxial ophthalmoscope, operating otoscope OR pneumatic otoscope with non-rechargeable or lithium ion 60-minute power handle and hard storage case

This Diagnostic Zip up Set makes portable diagnostics easy. The compact case has sections to keep all parts organised, whilst being quick and easy to access. Carry into any room, use outside of your practice, and share with co-workers whilst ensuring you never misplace a part again.

#### Choose from:

Operating - Non-Rechargeable Set

Pneumatic - Lithium Ion - Rechargeable Set

Operating - Lithium Ion - Rechargeable Set

**Key Features:**

**Coaxial Ophthalmoscope:**

- 6 apertures for all general and specialist use
- 8 lenses -25 to +40 diopters for better resolution
- Rubber brow rest prevents scratching of eyeglasses
- Illuminated lens dial clearly identifies settings
- Coaxial optical system for shadow-free spot and easier entry into undilated pupil
- Crossed linear polarizing filter to reduce corneal reflection
- Red-free filter for use with any aperture
- 18 possible aperture combinations for greater versatility
- Sealed optics for better dust protection

**Operating Otoloscope:**

- Open system for convenience during procedures
- Rotatable lens and speculum for ease of use
- Fibre optic illumination for cool light

**Pneumatic Otoloscope:**

- Rotating lens with ultra-seal for pneumatic
- Fibre optic illumination for cool light

**Non rechargeable handle for C cells:**

- Sure grip, smooth, nonabrasive knurled finish ensures a no-slip grip
- Chrome-plated brass handles with durable, long-life construction
- Safety switch to prevent accidental draining of battery
- Metal head connectors for durability

Each diagnostic set also includes 3 x polypropylene specula and (4,5 & 7mm) and zip case.

**ADDITIONAL INFORMATION**

<b>Animal</b>	Small Animal
<b>SKU</b>	WA962-case-set
<b>Overall Product Dimensions (mm)</b>	11.7cmW x 19.7cmH x 5.5cmD
<b>Colour/Style</b>	Black
<b>Brand</b>	Welch Allyn
<b>Catagory Sub Heading</b>	Save 20% (Price includes saving)
<b>Digital</b>	No

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights