

HANDI+ OXYGEN ANALYSER

129-555



PRODUCT DESCRIPTION

So you have an oxygen therapy cage, and the oxygen concentration you provide is?

If you can't answer this question, it's worth considering this product. Easily fitted to an O₂ therapy door the Handi+ displays the concentration level on an LCD display. It can also be used in line with your anaesthetic machine, ideal for verifying output from O₂ concentrators.

It is simple to use, with one-touch calibration, a new over-mold, calibration reminder, and auto-off, combined with a 2-year warranty to make it the best of its class.

Features

- New and improved one-touch calibration button
- Durable over-mold for extra protection
- Small, convenient, and portable
- Auto-off
- Lanyard Attachment
- Protective Wrap

ADDITIONAL INFORMATION

Animal	Mixed
SKU	129-555
Colour/Style	Blue
Monitor Parameters	Oxygen

Screen Type	Digital Display
Power Source	Battery
Optional Modules	No
SD Card	No

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.
Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!
Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.
Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights