

SILICON PROTECTIVE CASE FOR BP DOPPLER

DOP-400

**PRODUCT DESCRIPTION****Protect your Burtons BP Doppler with this Silicon Protective Case**

- Designed for the Burtons BP Doppler
- Non - slip silicon
- Flexible wrap around design

ADDITIONAL INFORMATION

Animal	Small Animal
SKU	DOP-400
Colour/Style	Clear
Brand	Burtons
Type	Other
Monitor Fit	Vet BP Doppler

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights