

SILICONE PROTECTIVE CASE FOR VP DOPPLER

DOP-400



PRODUCT DESCRIPTION

Protect your Burtons BP Doppler with this Silicon Protective Case

- Designed for the Burtons BP Doppler
- Non - slip silicon
- Flexibile wrap around design

Silicone Protective Case For VP Doppler. Change when broken. Compatible with VetBP monitors.

ADDITIONAL INFORMATION

Animal	N/A
SKU	DOP-400
Accessory Parameter	NIBP
Monitor Compatibility	Burtons PM-60, Mindray ePM12, Mindray uMEC
Change Frequency	When Broken

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights