

BURTONS CO² MONITOR EXHAUST KIT

V1175

**PRODUCT DESCRIPTION****Enable sample gases to be evacuated into your scavenge system**

- CO₂ Monitor Exhaust Kit
- Complete with Scavenging tube (Adaptor V1175-003) and Scavenge Connector 22mm M/F with port (V1175-003)
- For monitors: V8200, V8401, V8402, V90041, V90043 V9400 SurgiVet Advisor Monitor & All IMEC's

Change frequency when broken. Compatible with M880B, Edan X12, Mindray uMEC, Mindray ePM12 & SurgiVet Advisor monitors.

ADDITIONAL INFORMATION

Animal	N/A
SKU	V1175
Brand	Burtons
Accessory Parameter	CO ²
Monitor Compatibility	BLT M880B, Edan X12, Mindray ePM12, Mindray uMEC, SurgiVet Advisor
Change Frequency	When Broken

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights