

## BURTONS HANDHELD PULSE OXIMETER

129-810



### PRODUCT DESCRIPTION

**Intended for use on dogs, cats or similar size patients the durable, compact and lightweight Pulse Oximeter is ideal for continuous monitoring of oxygen saturation (SpO<sub>2</sub>)**

The compact handheld capnography, with its bright colour display and 100 patient storage, it can be used in any environment where continuous, non-invasive monitoring of these parameters is required.

**Limited stock remaining!**

#### Key Features

- Bright 2.4" TFT LCD Colour Display
- Suitable for small and large animals
- SpO<sub>2</sub> with waveforms
- Data storage for up to 100 patients over 72 hours.
- Audible and visual physiological and technical alarms
- Trend graph and trend table review
- Alarm event records review
- Internal battery
- External power supply

#### ADDITIONAL INFORMATION

<b>Animal</b>	Small Animal
<b>SKU</b>	129-810
<b>Colour/Style</b>	White
<b>Brand</b>	Burtons
<b>Monitor Parameters</b>	SpO <sub>2</sub>
<b>Screen Type</b>	Full Colour
<b>Power Source</b>	Battery, Rechargeable
<b>Optional Modules</b>	No
<b>SD Card</b>	No

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. **Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! **Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. **Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights